



**Material Safety Data Sheet**  
**Carbon Standards, Organic, COD Standards,**  
**Potassium Acid Phthalate, 0.05 N**

**Section 1 - Chemical Product and Company Identification**

**MSDS Name:**

Carbon Standards, Organic, COD Standards, Potassium Acid Phthalate, 0.05 N

**Catalog Numbers:**

LC12910, LC13240, LC13245, LC18700

**Synonyms:**

Potassium hydrogen phthalate solution

**Section 2 - Composition, Information on Ingredients**

CAS#	Chemical Name:	Percent
7732-18-5	Water	balance
877-24-7	Potassium hydrogen phthalate	<1
7664-38-2	Phosphoric acid	<1

**Section 3 - Hazards Identification**

**Emergency Overview**

**Appearance:** Clear, colorless solution

**Caution!** May cause irritation to eyes, skin, respiratory and gastrointestinal tracts. May be harmful if swallowed.

**Target Organs:** Eyes, skin, respiratory and gastrointestinal tracts

**Potential Health Effects**

**Eye:**

May cause eye irritation.

**Skin:**

May cause skin irritation.

**Ingestion:**

Causes gastrointestinal irritation with nausea, vomiting and diarrhea.



**Material Safety Data Sheet**  
**Carbon Standards, Organic, COD Standards,**  
**Potassium Acid Phthalate, 0.05 N**

**Inhalation:**

May cause respiratory tract irritation.

**Chronic:**

May cause dermatitis and conjunctivitis. Repeated inhalation may cause lung damage.

**Section 4 - First Aid Measures**

**Eyes:**

Flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower lids until no evidence of chemical remains. Get medical aid at once.

**Skin:**

Flush skin with plenty of soap and water for at least 15 minutes while removing contaminated clothing and shoes. Get medical aid.

**Ingestion:**

Give conscious victim 2-4 cupfuls of milk or water. Never give anything by mouth to an unconscious person. Get medical aid at once.

**Inhalation:**

Move victim to fresh air immediately. If breathing is difficult, administer oxygen. Give artificial respiration if necessary. Get medical aid at once.

**Notes to Physician:**

Treat symptomatically and supportively.

**Section 5 - Fire Fighting Measures**

**General Information:**

As in any fire, wear a self-contained breathing apparatus in pressure-demand, MSHA/NIOSH (approved or equivalent), and full protective gear. During a fire, irritating and highly toxic gases may be generated by thermal decomposition or combustion.

**Extinguishing Media:**

For small fires, use dry chemical, carbon dioxide, water spray or alcohol-resistant foam.

**Autoignition Temperature:**

Not available

**Flash Point:**

Not applicable

**NFPA Rating:**

Health- 1, Flammability- 0, Instability- 0

**Explosion Limits:**

Lower: n/a      Upper: n/a

**Section 6 - Accidental Release Measures**

**General Information:**

Use proper personal protective equipment as indicated in Section 8.

**Spills/Leaks:**

Absorb spill with inert material such as sand, vermiculite, or diatomaceous earth, and transfer to a suitable container labeled for later disposal.



**Material Safety Data Sheet**  
**Carbon Standards, Organic, COD Standards,**  
**Potassium Acid Phthalate, 0.05 N**

**Section 7 - Handling and Storage**

**Handling:**

Wash thoroughly after handling. Do not ingest or inhale. Do not get in eyes, on skin, or on clothing.

**Storage:**

Store tightly capped in a cool, dry, well-ventilated area away from incompatible materials.

**Section 8 - Exposure Controls, Personal Protection**

**Engineering Controls:**

Facilities using or storing this material should be equipped with an eyewash and safety shower. Provide local exhaust or general dilution ventilation to keep airborne levels below the permissible exposure limits.

**Exposure Limits:**

<b>Chemical Name:</b>	<b>ACGIH</b>	<b>NIOSH</b>	<b>OSHA</b>
<b>Water</b>	None listed	None listed	None listed
<b>Potassium hydrogen phthalate</b>	None listed	None listed	None listed
<b>Phosphoric acid</b>	1 mg/m3 TWA 3 mg/m3 STEL	1 mg/m3 TWA 1000 mg/m3 IDLH	1 mg/m3 TWA

**OSHA Vacated PELs:**

Phosphoric acid: 1 mg/m3 TWA

**Personal Protective Equipment**

**Eyes:**

Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133. Do not wear contact lenses when working with chemicals.

**Skin:**

Wear appropriate protective gloves to prevent skin exposure.

**Clothing:**

Wear appropriate protective clothing to prevent skin exposure.

**Respirators:**

Follow the OSHA respirator regulations found in 29 CFR 1910.134. Always use a NIOSH-approved respirator when necessary.

**Section 9 - Physical and Chemical Properties**

**Physical State:** Liquid  
**Color:** Colorless  
**Odor:** Odorless  
**pH:** Acidic



## Material Safety Data Sheet Carbon Standards, Organic, COD Standards, Potassium Acid Phthalate, 0.05 N

<b>Vapor Pressure:</b>	Not available
<b>Vapor Density:</b>	Not available
<b>Evaporation Rate:</b>	Not available
<b>Viscosity:</b>	Not available
<b>Boiling Point:</b>	Not available
<b>Freezing/Melting Point:</b>	Not available
<b>Decomposition Temperature:</b>	Not available
<b>Solubility in water:</b>	Soluble
<b>Specific Gravity/Density:</b>	1.0
<b>Molecular Formula:</b>	Not available
<b>Molecular Weight:</b>	Not available

### Section 10 - Stability and Reactivity

**Chemical Stability:**

Stable under normal temperatures and pressures.

**Conditions to Avoid:**

Incompatible materials, excess heat.

**Incompatibilities with Other Materials:**

Strong oxidizing agents, bases, alcohols, amines, metals, aldehydes, reducing agents, organic peroxides, ketones, halogenated hydrocarbons.

**Hazardous Decomposition Products:**

Carbon monoxide, carbon dioxide, potassium oxides, potassium fume, hydrogen gas, oxides of phosphorus, phosphine

**Hazardous Polymerization:**

Has not been reported.

### Section 11 - Toxicological Information

**RTECS:**

CAS# 7732-18-5: ZC0110000

CAS# 877-24-7: CZ4326000

CAS# 7664-38-2: TB6300000

**LD50/LC50:**

CAS# 7732-18-5:

Oral, rat: LD50 = >90 mL/kg

CAS# 877-24-7:

Oral, rat: LD50 = >3200 mg/kg

Skin, guinea pig: LD50 = >1 g/kg

CAS# 7664-38-2:

Inhalation, rat: LC50 = 25.5 mg/m<sup>3</sup>

Oral, rat: LD50 = 1530 mg/kg

Skin, rabbit: LD50 = 2740 mg/kg

**Carcinogenicity:**

CAS# 7732-18-5: Not listed as a carcinogen by ACGIH, IARC, NTP, OSHA, or California Proposition 65.



## Material Safety Data Sheet Carbon Standards, Organic, COD Standards, Potassium Acid Phthalate, 0.05 N

CAS# 877-24-7: Not listed as a carcinogen by ACGIH, IARC, NTP, OSHA, or California Proposition 65.

CAS# 7664-38-2: Not listed as a carcinogen by ACGIH, IARC, NTP, OSHA, or California Proposition 65.

**Epidemiology:**

See actual entry in RTECS for complete information.

**Teratogenicity:**

Teratogenic effects have occurred in laboratory animals.

**Reproductive:**

Experiments have shown reproductive toxicity effects on laboratory animals.

**Mutagenicity:**

Mutagenic effects have occurred in microorganisms.

**Neurotoxicity:**

No information found.

### Section 12 - Ecological Information

No information found.

### Section 13 - Disposal Considerations

Dispose of in accordance with Federal, State, and local regulations.

### Section 14 - Transport Information

**US DOT**

**Shipping Name:** Not regulated

**Hazard Class:**

**UN Number:**

**Packing Group:**

### Section 15 - Regulatory Information

**US Federal**

**TSCA:**

CAS# 7732-18-5 is listed on the TSCA Inventory.

CAS# 877-24-7 is listed on the TSCA Inventory.

CAS# 7664-38-2 is listed on the TSCA Inventory.

**SARA Reportable Quantities (RQ):**

CAS# 7664-38-2: final RQ = 5000 pounds (2270 kg)

**CERCLA/SARA Section 313:**

Not reportable under Section 313

**OSHA - Highly Hazardous:**

Not considered highly hazardous by OSHA



**Material Safety Data Sheet**  
**Carbon Standards, Organic, COD Standards,**  
**Potassium Acid Phthalate, 0.05 N**

**US State**

**State Right to Know:**

Phosphoric acid can be found on the following state Right-to-Know lists: California, New Jersey, Florida, Pennsylvania, Minnesota, Massachusetts.

**California Regulations:**

None

**European/International Regulations**

**Canadian DSL/NDSL:**

CAS# 7732-18-5 is listed on Canada's DSL List.

CAS# 877-24-7 is listed on Canada's DSL List.

CAS# 7664-38-2 is listed on Canada's DSL List.

**Canada Ingredient Disclosure List:**

CAS# 7732-18-5 is not listed on Canada's Ingredient Disclosure List.

CAS# 877-24-7 is not listed on Canada's Ingredient Disclosure List.

CAS# 7664-18-5 is listed on Canada's Ingredient Disclosure List.

**Section 16 - Other Information**

MSDS Creation Date: October 8, 2010

Revision Date: October 11, 2010

*Information in this MSDS is from available published sources and is believed to be accurate. No warranty, express or implied, is made.*